



F-900 Controller Software and Calibration Instructions

Pre-Calibration Preparation

The following is a listing of the total standard gases/materials needed to calibrate ALL sensors available in the F-900 (Not all F-900's come equipped with all sensors, see graph to verify which calibration gases to purchase for the sensors in your instrument) **Every F-900** is equipped with C₂H₄ ppm sensor and C₂H₄ ppb sensor , **optional sensors** are: CO₂ PCT sensor , CO₂ ppm sensor, O₂ sensor:

- Ethylene calibration gas for C₂H₄ ppm sensor (20 ppm recommended)
- Ethylene verification gas for C₂H₄ ppm sensor (10 ppm recommended)
- Ethylene calibration gas for C₂H₄ ppb sensor (1.5 ppm recommended)
- Ethylene verification gas for C₂H₄ ppb sensor (1 ppm recommended)
- CO₂ calibration gas for CO₂ PCT sensor (16% recommended)
- CO₂ verification gas for CO₂ PCT sensor (5% recommended)
- CO₂ calibration gas for CO₂ ppm sensor (2,000 ppm recommended)
- CO₂ verification gas for CO₂ ppm sensor (1,000 ppm recommended)
- O₂ calibration gas (21% recommended)* can use ambient air instead
- 100% Nitrogen (N₂) gas for CO₂ PCT sensor and O₂ sensor
- Potassium Permanganate (KMnO₄) granules
- Soda Lime (NaOH, KOH, Ca(OH)₂) granules
- Flow Meter Capable of Measuring 300 mL/min and lower
- Tubing with low outgassing properties (Viton)

*The listed gases can be replaced with similar concentrations as long as the gas used is of a standard, consistent value. The following is a chart displaying the use of each gas to the appropriate sensor:

Standard Gas Needed for Each Sensor:

	20 ppm (C ₂ H ₄) gas	10 ppm (C ₂ H ₄) gas	1.5 ppm (C ₂ H ₄) gas	1.0 ppm (C ₂ H ₄) gas	16% (CO ₂) gas	5% (CO ₂) gas	2,000 ppm (CO ₂) gas	1,000 ppm (CO ₂) gas	21% (O ₂) gas **	100% (N ₂) gas	Flow Meter (≤300 mL/min)	Soda Lime (NaOH-KOH...)	Potassium Permanganate (KMnO ₄)
C ₂ H ₄ ppm sensor	✓	✓									✓		✓
C ₂ H ₄ ppb sensor			✓	✓							✓		✓
CO ₂ PCT sensor					✓	✓				✓	✓		
CO ₂ ppm sensor							✓	✓		✓	✓	✓	
O ₂ sensor									✓	✓	✓		

** Can use ambient air instead of purchasing 21% (O₂) gas

Calibration Schedule for F-900:

Calibration	C ₂ H ₄ ppm sensor	C ₂ H ₄ ppb sensor	CO ₂ PCT sensor	CO ₂ ppm sensor	O ₂ sensor
Set Zero	Daily	Daily	6 months	Weekly	6 months
Set Span	6 months	Weekly	6 months	6 months	6 months

Controller Software

*not required for calibrations involving Potassium Permanganate (KMnO₄):

1. If you have not already done so, download the package setup software from <http://www.felixinstruments.com/support/f-900-support/software> and install it on your computer. This package includes the latest firmware updates for the F-900, as well as the controller software needed for real time monitoring and calibration of the F-900.
2. Turn on the F-900.
3. Connect the F-900 to computer via a USB cable.
4. On the F-900, from the main menu, go to Setup > Measure and change the “Mode” setting to “Monitor”. Return to the main menu by pressing the left arrow.
5. Open the controller software.
6. In the top left of the window, select the F-900 unit listed.
7. Click the “Calibration” tab and select the “900” sub-tab.

Ethylene (C₂H₄) PPB Sensor Calibration:


Setting the Zeroes (Offset) with Potassium Permanganate (KMnO₄) *daily

1. Place fresh potassium permanganate (KMnO₄) in Chamber Out (Solids Chamber).
2. Power on the instrument and allow adequate warm-up time (3 min or more).
3. Press the Right arrow when Measure is highlighted to begin a measurement.
4. Let the sensor automatically “Environment Stabilize”
5. Once the measurement begins, press the right arrow until you access the “settings” menu.
6. Scroll down to “Set Zero”
7. Press the Right arrow to highlight “C₂H₄”
8. Press Enter.
9. A message appears asking “zero selected sensors?”
 - a. Press enter to confirm.
10. A message appears asking to “Place KMnO₄ in CH_Out.”
 - a. If KMnO₄ is in Chamber Out, press enter to confirm
11. The display will switch to Monitor mode and “correcting offset” is shown at the bottom.
 - a. The settings will automatically change to turn on Chamber Out, turn off Chamber In and turn on the Close loop.

12. The instrument will wait until no more than a 40 ppb change is detected for 10 minutes before setting the zero calibration point.
13. The instrument will make a “beep” sound twice to indicate that the offset correction is complete. The entire set zero process takes about 25 minutes.
14. The measurement will continue with the original settings for the conditioning chambers and close loop on/off.

Ethylene (C₂H₄) PPB Sensor Set Span *Weekly

1. Connect the high resolution flowmeter, ethylene calibration gas for the ppb sensor, and F-900 intake ports together using a T-junction of tubing.
2. On the F-900, begin Measure mode.
3. Set the flowrate on the C₂H₄ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.
4. Let the F-900 read the gas for five minutes.
5. In the controller software, in the C₂H₄ field, select PPB and set the concentration to the value of the calibration gas being used. Then hit the “Set Span” button.



6. The concentration on the F-900 should now read the same as the calibration gas 7% ± 25 PPB.

**An alternate method for C₂H₄ PPB span calibration is available with purchase of the F-900 calibration kit. For instructions on this method, see pg. 44 of the F-900 user manual entitled “C₂H₄ PPB Re-Calibration (Manual)”.*

Ethylene (C₂H₄) PPB Sensor Verification *After Set Span

Perform this process immediately after calibrating the ethylene ppb sensor on the F-900.

1. Connect the high resolution flowmeter, ethylene verification gas for the ppb sensor, and F-900 intake ports together using a T-junction of tubing.
2. On the 900 begin measure mode.
3. Set the flowrate on the C₂H₄ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the F-900 again. If the error persists, contact tech support.

4. Let the F-900 read the gas for five minutes.
 5. The F-900 should read the same value as the verification gas within $7\% \pm 25$ PPB . If it does not, repeat the “Ethylene Sensor Calibration” section above.
-

Ethylene (C₂H₄) PPM Sensor Calibration


Setting the Zeroes (Offset) with Potassium Permanganate (KMnO₄) *daily

1. Place fresh potassium permanganate (KMnO₄) in Chamber Out.
2. Power on the instrument and allow adequate warm-up time (3 min).
3. Press the Right arrow when Measure is highlighted to begin a measurement.
4. Let the sensor automatically “Environment Stabilize”
5. Once the measurement begins, press the right arrow until you access the “settings” menu.
6. Scroll down to “Set Zero”.
7. Press the Right arrow to highlight “C₂H₄”.
8. Press Enter.
9. A message appears asking “zero selected sensors?”
 - a. Press enter to confirm.
10. A message appears asking to “Place KMnO₄ in CH_Out.”
 - a. If KMnO₄ is in Chamber Out, press enter to confirm
11. The display will switch to Monitor mode and “correcting offset” is shown at the bottom.
 - a. The settings will automatically change to turn on Chamber Out, turn off Chamber In and turn on the Close loop.
12. The instrument will wait until no more than a 40 ppb change is detected for 10 minutes before setting the zero calibration point.
13. The instrument will make a “beep” sound twice to indicate that the offset correction is complete. The entire set zero process takes about 25 minutes.
14. The measurement will continue with the original settings for the conditioning chambers and close loop on/off.

Ethylene (C₂H₄) PPM Sensor Set Span *6 months

1. Connect the high resolution flowmeter, ethylene calibration gas for the ppm sensor, and F-900 intake ports together using a T-junction of tubing.
2. On the F-900, begin measure mode.
3. Set the flowrate on the C₂H₄ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.

- Let the F-900 read the gas for five minutes.
- In the controller software, in the C₂H₄ field, select PPM and set the concentration to the value of the calibration gas being used. Then hit the “Set Span” button.



- The concentration on the 900 should now read the same as the calibration gas 5% ± 0.5 PPM .

Ethylene (C₂H₄) PPM Sensor Verification *After Set Span

Perform this process immediately after calibrating the ethylene ppm sensor on the F-900.

- Connect the high resolution flowmeter, Ethylene verification gas for the ppm sensor, and F-900 intake ports together using a T-junction of tubing.
- On the F-900 begin measure mode.
- Set the flowrate on the C₂H₄ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the F-900 again. If the error persists, contact tech support.
- Let the F-900 read the gas for five minutes.
- The F-900 should read the same value as the verification gas within 5% ± 0.5 PPM. If it does not, repeat the “Ethylene Sensor Calibration” section above.


Carbon Dioxide (CO₂) PCT Sensor Calibration

*CO₂ PCT Sensor is optional

Setting the Zeroes (Offset) with Nitrogen Gas (N₂) *6 months

- Connect the high resolution flow meter, N₂ gas, and F-900 intake ports together using a T-junction of tubing.
- On the F-900, begin measure mode.
- Set the flowrate on the N₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the F-900 again. If the error persists, contact tech support.
- Let the F-900 read the gas for five minutes.


- In the controller software, select CO₂ PCT in the “Zero sensors” field, then click the “Set Zero” button.



- The F-900 should now read 0 for the CO₂ PCT sensor.
- Set Span should take place directly following the Set Zero process for the CO₂ PCT sensor!

Carbon Dioxide (CO₂) PCT Sensor Set Span *6 months

- Set Span should take place directly following the Set Zero process for the CO₂ PCT sensor.
- Connect the high resolution flow meter, CO₂ calibration gas for the PCT sensor, and 900 intake ports together using a T-junction of tubing.
- On the F-900, begin measure mode.
- Set the flowrate on the CO₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.
- Let the F-900 read the gas for five minutes.
- In the controller software, in the CO₂ field, select PCT and set the concentration to the value of the calibration gas being used. Then hit the “Set Span” button.



- The concentration on the 900 should now read the same as the calibration gas \pm 5%.

Carbon Dioxide (CO₂) PCT Sensor Verification *After Set Span

Perform this process immediately after calibrating the carbon dioxide sensor on the F-900

- Connect the high resolution flow meter, CO₂ PCT sensor verification gas, and F-900 intake ports together using a T-junction of tubing.
- On the F-900, begin measure mode.
- Set the flowrate on the CO₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.
- Let the F-900 read the gas for five minutes.


5. The F-900 should read the same value as the verification gas \pm 5%. If it does not, repeat the “Carbon Dioxide PCT Sensor Calibration” section above.

Carbon Dioxide (CO₂) PPM Sensor Calibration

*CO₂ ppm Sensor is optional

Setting the Zeroes (Offset) with Soda Lime or Nitrogen Gas (N₂) *Weekly


1. Place granular Soda Lime in the solids chamber of the F-900, or if using nitrogen gas, connect the high resolution flow meter, N₂ gas, and F-900 intake ports together using a T-junction of tubing.
2. On the F-900, begin measure mode.
3. If using Soda Lime, run the instrument on closed loop, otherwise, set the flowrate on the N₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.
4. Let the F-900 measure for five minutes.
5. In the controller software, select CO₂ ppm in the “Zero sensors” field, then click the “Set Zero” button (it will display “Use N₂” when utilizing both Soda Lime or N₂ gas).



6. The 900 should now read 0 for the CO₂ ppm sensor.

Carbon Dioxide (CO₂) PPM Sensor Set Span *6 months

1. Connect the high resolution flow meter, CO₂ calibration gas for the ppm sensor, and F-900 intake ports together using a T-junction of tubing.
2. On the F-900, begin measure mode.
3. Set the flowrate on the CO₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the F-900 again. If the error persists, contact tech support.
4. Let the F-900 read the gas for five minutes.
5. In the controller software, in the CO₂ field, select ppm and set the concentration to the value of the calibration gas being used. Then hit the “Set Span” button.



6. The concentration on the F-900 should now read the same as the calibration gas $\pm 5\%$.

Carbon Dioxide (CO₂) PPM Sensor Verification *After Set Span

Perform this process immediately after calibrating the carbon dioxide PPM sensor on the F-900.

1. Connect the high resolution flow meter, CO₂ ppm sensor verification gas, and F-900 intake ports together using a T-junction of tubing.
2. On the F-900, begin measure mode.
3. Set the flowrate on the CO₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the F-900 again. If the error persists, contact tech support.
4. Let the F-900 read the gas for five minutes.
5. The F-900 should read the same value as the verification gas $\pm 5\%$. If it does not, repeat the “Carbon Dioxide ppm Sensor Calibration” section above.

Oxygen (O₂) Sensor Calibration

*Oxygen sensor is optional

Setting the Zeroes (Offset) with Nitrogen Gas (N₂) *6 months

1. Connect the high resolution flow meter, N₂ gas, and F-900 intake ports together using a T-junction of tubing.
7. On the F-900, begin measure mode and measure for 5 minutes
8. Set the flowrate on the N₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.
9. Let the F-900 read the gas for five minutes.
10. On the F-900, navigate to Setup > Mode > USB > Controller, to enable connection to the controller software
11. Connect the F-900 to a computer via USB and launch the controller software
1. Select the appropriate F-900 serial number in the upper left of the controller software screen

V0.7 Customer Controller Software Instructions


- In the controller software, select O₂ in the “Zero sensors” field, then click the “Set Zero” button.



- The F-900 should now read 0 for the O₂ sensor.

Oxygen (O₂) Sensor Set Span *6 months

- Connect the high resolution flow meter, O₂ calibration gas, and F-900 intake ports together using a T-junction of tubing.
- On the F-900 begin measure mode.
- Set the flowrate on the O₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the 900 again. If the error persists, contact tech support.
- Let the F-900 read the gas for five minutes.
- In the controller software, in the O₂ field, set the concentration to the value of the calibration gas being used. Then hit the “Set Span” button.



- The concentration on the F-900 should now read the same as the calibration gas \pm 5%.

Oxygen (O₂) Sensor Verification *After Set Span

Perform this process immediately after calibrating the oxygen sensor on the F-900.

- Connect the high resolution flow meter, O₂ verification gas, and F-900 intake ports together using a T-junction of tubing.
- On the F-900 begin measure mode.
- Set the flowrate on the O₂ gas tank so that the flow meter reads between 70 and 200 mL/min. If the F-900 gives a “Flow Blocked” error, then increase the flow on the meter by 10 mL/min and begin measure mode on the F-900 again. If the error persists, contact tech support.
- Let the F-900 read the gas for five minutes.
- The F-900 should read the same value as the verification gas \pm 5%. If it does not, repeat the “Oxygen Sensor Calibration” section above.

Monitor Mode Verification

The following procedure is to check whether the sensors were calibrated to read within spec of the actual gas value.

It is highly recommended that this step be performed after one full day has passed since the calibration was performed. If this step is not performed, accuracy of the calibration cannot be verified.

Repeat the verification sections from above for each gas and record the results here after three minutes (at least) of measure time. If any of the values are out of specification, check your procedure and recalibrate the sensor again.

Gas	900 Results	Spec: $\pm 5\%$
Air (Air is ~20.9% O₂)		<input type="checkbox"/> Yes <input type="checkbox"/> No
CO₂ PCT Verification Gas		<input type="checkbox"/> Yes <input type="checkbox"/> No
CO₂ PPM Verification Gas		<input type="checkbox"/> Yes <input type="checkbox"/> No
O₂ Verification Gas		<input type="checkbox"/> Yes <input type="checkbox"/> No

Gas	900 Results	Spec: 7%, ± 25 PPB
C₂H₄ PPB Verification Gas		<input type="checkbox"/> Yes <input type="checkbox"/> No

Gas	900 Results	Spec: 5%, ± 0.5 PPM
C₂H₄ PPM Verification Gas		<input type="checkbox"/> Yes <input type="checkbox"/> No

For information on sourcing known gases required for calibration, please refer to the following website to inquire about your region: <https://www.airliquide.com/group/where-we-operate>