

F-900

■ Accessories



F-900IK ■ Injection Port Kit

The Injection Port Kit allows the F-900 to operate like a Gas Chromatograph and analyze low concentrations in small-volume samples.



F-900CK ■ Calibration Kit

Includes materials to perform calibrations on the F-900 Portable Ethylene Analyzer PPB ethylene sensor. We recommend that span calibration of this high-sensitivity sensor be performed weekly to ensure that the accuracy of the device is maintained.



F-900RK ■ Research Kit

The F-900 Research Kit provides a comprehensive way to measure directly from fruits or from static samples in any environment.

The F-900RK includes the following (F-900 Ethylene Analyzer sold separately):

- ▶ 2 liter chamber with solid bottom plate for closed air-tight measurement, and bottom plate with seal for on-plant measurement
- ▶ Chamber wand with flexible tubing assembly
- ▶ Syringe
- ▶ Sampling bag
- ▶ Silica gel
- ▶ Injection Port Kit for small volume sampling
- ▶ Hard shell instrument carrying case