

Accessories



LC-1
■ Square Leaf Chamber

For open-system measurements of trees, shrubs and herbs with small, broad leaves.
25 mm x 25 mm



LC-2
■ Wide Rectangular Leaf Chamber

For open-system measurements of trees, shrubs and herbs with small, regular leaves.
55 mm x 20 mm



LC-3
■ Narrow Rectangular Leaf Chamber

For open-system measurements of grasses and grass-like leaves.
65 mm x 10 mm



LC-4
■ Small Cylindrical Leaf Chamber

For open-system measurements of succulents and small-needled conifers.
25 mm x 90 mm



LC-5
■ Large Cylindrical Leaf Chamber

For open-system measurements of large-needled conifers.
50 mm x 70 mm



LC-7
■ 1/4 Liter Leaf Chamber

For closed-system measurements of long, narrow leaves.
104 mm x 33 mm x 73 mm



LC-8
■ 1/2 Liter Leaf Chamber

For closed-system measurements of medium-sized leaves.
89 mm x 66 mm x 86 mm



LC-9
■ 1 Liter Leaf Chamber

For closed-system measurements of large leaves.
112 mm x 91 mm x 99 mm



LC-10
■ 4 Liter Leaf Chamber

For closed-system measurements of very large leaves.
180 mm x 130 mm x 170 mm



LC-11
■ Cactus Leaf Chamber

For measuring the leaves of Cacti with the CI-340 Handheld Photosynthesis System.
Window size: 14.6cm²



CI-301SR
■ Soil Respiration Chamber

For measurements of CO₂ flux from soils.



CI-301CC
■ Canopy Chamber Attachment

Interface between chamber and instrument which allows user to design and build their own custom chamber.



F-9002L
■ Sample Chamber

For measuring small fruits. Includes a solid bottom plate for closed air-tight measurement.



F-9008L
■ Sample Chamber

Air-tight chamber for measuring large fruits.



F-900CW
■ Attachment Wand

Interface for the F-900 and any sample chamber. Also included in F-900RK Research Kit.